


Revo-S™ SC1 SC2 SU

Protocolle Protocol Protokoll



Vitesse 250 à 400 tr/min - Couple maxi 0.8 N.cm Speed 250 to 400 rpm - Maxi Torque 0.8 N.cm
Motorgeschwindigkeit: 400 U/min - Max. Drehmoment: 0.8 N.cm



www.micro-mega.com



Your Endo Specialist™

Revo-S™ AS30 AS35 AS40

Protocole Protocol Protokoll

AS30

N°30 6%
L 21 - 25 -
29 mm



LT-WL-AL

AS35

N°35 6%
L 21 - 25 -
29 mm



LT-WL-AL

AS40

N°40 6%
L 21 - 25 -
29 mm



LT-WL-AL

Vitesse 250 à 400 tr/min - Couple maxi 0.8 N.cm Speed 250 to 400 rpm - Maxi Torque 0.8 N.cm
Motorgeschwindigkeit: 400 U/min - Max. Drehmoment: 0.8 N.cm



Your Endo Specialist™

www.micro-mega.com

